

Anastomosis Clip


Mizuho Anastomosis Clips are dependable and made of Elgiloy for superior durability and spring elasticity, like their predecessor, the Sugita Aneurysm Clips.

Both high and low pressure clips are available with a wide variety of sizes and styles to accommodate various surgical approaches.

Closing force is measured 1.0mm from the distal tip of the blade when opened 1.0mm wide by computerized electronic gauge measuring equipment.



Anastomosis Clip



Prod. No.	Product Description	b	L	Closing force
07-096-11	Anastomosis Clip Straight 2mm LP	2mm	6mm	20-40g 0.20-0.39N
07-096-12	Anastomosis Clip Straight 3mm LP	3mm	7mm	
07-096-13	Anastomosis Clip Straight 4mm LP	4mm	8mm	
07-096-14	Anastomosis Clip Straight 5mm LP	5mm	9mm	
07-096-15	Anastomosis Clip Bent 3mm LP	3mm	7mm	
07-096-16	Anastomosis Clip Bent 4mm LP	4mm	8mm	
07-096-17	Anastomosis Clip Sideward Bent LP	3mm	7.9mm	50-70g 0.49-0.69N
07-096-21	Anastomosis Clip Straight 2mm HP	2mm	6mm	
07-096-22	Anastomosis Clip Straight 3mm HP	3mm	7mm	
07-096-23	Anastomosis Clip Straight 4mm HP	4mm	8mm	
07-096-24	Anastomosis Clip Straight 5mm HP	5mm	9mm	
07-096-25	Anastomosis Clip Bent 3mm HP	3mm	7mm	
07-096-26	Anastomosis Clip Bent 4mm HP	4mm	8mm	
07-096-27	Anastomosis Clip Sideward Bent HP	3mm	7.9mm	


b = blade length L = total length

Dual Needle Holder

Two-in-One Instrument

Dual Needle Holders also act as micro forceps, so no need to switch between the two instruments making surgical procedures smoother.


Prod. No.	Product Description	Length
07-861-04	Dual Needle Holder Straight	13cm
07-861-05	Dual Needle Holder Curved	13cm
07-861-14	Dual Needle Holder Straight (For Deep Area)	18cm
07-861-15	Dual Needle Holder Curved (For Deep Area)	18cm



Anastomosis Clip Applier

Two Types of Appliers

- Standard, Bayonet applier in 17.5cm length
- Straight, Short applier in 14.7cm length



Prod. No.	Product Description	Length
07-097-00	Anastomosis Clip Applier	17.5cm
07-097-01	Anastomosis Clip Applier Straight / Short	14.7cm

* No CE Marking

Mizuho Medical Co., Ltd.

Hongo Shintoku Building 7F, 3-38-1
Hongo, Bunkyo-ku, Tokyo 113-0033, JAPAN

Phone +81-3-3815-7101
F a x +81-3-3815-1705
U R L www.mizuhomedical.co.jp
E-mail sales-dept@mizuhomedical.co.jp

Mizuho America Inc.

30057 Ahern Avenue Union City, CA 94587 Phone No.510 324 4500

Manufacturer

MIZUHO Corporation

Tokyo, JAPAN

Note: Published specifications are subject to change without prior notice.
Details and colors of the images might be different from ones of real products.

 Caution	<ul style="list-style-type: none"> ● Prior to product use, carefully read and comprehend the contents of the user manual and product packaging information.
	<ul style="list-style-type: none"> ● Inspect the product routinely and thoroughly, in addition to pre-and postoperative examinations.
	<ul style="list-style-type: none"> ● For after-sale support, contact your local sales representative.